

MOUTHGUARDS

Mouthguards should be worn during playing and training for any sport which could possibly contact the facial structures. The Australian Dental Association Inc (ADA) advises only using a professionally custom-made and custom-fitted mouthguard. Southern Cross Dental offers a wide range of colours and designs to choose from.

Colour Chart



White



Black



Yellow



Orange



Red



Pink



Lilac



Blue



Silver



Gold



Neon-Green



Neon-Yellow



Neon-Red

* Thickness : Transparent blank 2.0, 3.0, 4.0, 5.0 mm; Colour blank only 3.0 mm
*Photos for illustration purposes only. Colours may vary. Colours and labels attract extra charges.



YBY



WBW



RBR



BRY



BWR



RGB



Glitter-Yellow



Glitter-Orange



Glitter-Pink



Glitter-Purple



Glitter-Blue



Glitter-Green

Specifications

Junior

ALSO AVAILABLE IN QUEST



Thickness
3mm



Applications
Low impact sports
(for ages 5-12).

Examples include:
Speed skating, volleyball, soccer, little league sports.

Standard

ALSO AVAILABLE IN QUEST



Thickness
5mm



Applications
Non-contact sports (adults).

Examples include:
Intense cycling, rock climbing and weightlifting.

Professional

ALSO AVAILABLE IN QUEST



Thickness
6mm + 3mm
hard labial insert



Applications
Heavy contact sports where hard, pointed blows may be received by the jaw.

Examples include:
Hockey, rugby, boxing and wrestling.

QUEST RANGE

Mouthguard Colour Chart

Plain Colours



Red



Yellow



Green



Blue



Purple



Maroon



Baby Blue



Teal



Pink



Black



Brown



Navy Blue



Grey



Clear



White



Silver



Gold



Speed Sparks

Fluorescent Colours



Fluoro Green



Fluoro Yellow



Fluoro Orange



Fluoro Red

Multi-Coloured



Blue/Yellow



Blue/White



Red/White



Blue/Red



Red/Yellow



Black/Yellow



Black/White



Black/Red



Brown/Yellow



Purple/White



R/W/B



O/W/B



B/Y/M



B/Y/R



B/W/T



R/W/B

